

# Environmental & food safety products

## ENVIROMAX® (DRY) & ENVIROMAX PLUS® (PRE-MOISTENED)

Puritan	RELATED PRODUCTS				PRODUCT DESCRIPTION
	Copan	3M	Hardy Diagnostics	Medical Wire	
<b>88010-PF</b>	-	ENVSWB25D	88010PF	-	Pointed, dry foam tip, heat staked to polypropylene handle. Non-sterile.
<b>88010-PFB RND</b>	-	-	-	-	Rounded, dry foam tip, heat staked to polypropylene handle. Non-sterile.
<b>25-88050 PF</b>	-	ENVSWB25	2588050PF		Pointed, pre-moistened foam tip, heat staked to polypropylene handle. Sterile/individually wrapped.

## ESK® SAMPLING KIT: BUTTERFIELD'S SOLUTION

Puritan	RELATED PRODUCTS				PRODUCT DESCRIPTION
	Copan	3M	Hardy Diagnostics	Medical Wire	
<b>25-83004 PD BS</b>	-	-	2583004PDBS	MW781 w/5ml	Tube: Butterfield's Solution, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
<b>25-83004 PDB BS</b>	-	-	2583004PDBBS	-	
<b>25-83010 PD BS</b>	R4140 w10/ml	-	2583010PDBS	MW780 w/10ml	Tube: Butterfield's Solution, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
<b>25-83010 PDB BS</b>	-	-	2583010PDBBS	-	



**CALL** 800-321-2313 (USA & Canada) • 207-876-3311  
**EMAIL** sales@puritanmedproducts.com  
**VISIT** puritanmedproducts.com

PURITAN MEDICAL PRODUCTS CO.  
 GUILFORD, MAINE, 04443, USA  
 ISO 13485:2016 CE



**MADE IN THE USA**

## ESK® SAMPLING KIT: LETHEEN BROTH

Puritan	RELATED PRODUCTS				PRODUCT DESCRIPTION
	Copan	3M	Hardy Diagnostics	Medical Wire	
25-83004 PD LB	–	RS9604LET w/4ml	2583004PDLB	MW793 w/5ml	Tube: Letheen Broth, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83004 PDB LB	–	–	2583004PDBLB	–	
25-83010 PD LB	R4130 w/10ml	RS96010LET w/10ml	2583010PDLB	MW792 w/10ml	Tube: Letheen Broth, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83010 PDB LB	–	–	2583010PDBLB	–	

## ESK® SAMPLING KIT: NEUTRALIZING BUFFER

Puritan	RELATED PRODUCTS				PRODUCT DESCRIPTION
	Copan	3M	Hardy Diagnostics	Medical Wire	
25-83004 PD NB	–	RS9604NB w/4ml	2583004PDNB	MW784 w/5ml	Tube: Neutralizing Buffer, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83004 PDB NB	–	–	2583004PDBNB	–	
25-83010 PD NB	–	RS96010NB w/10ml	2583010PDNB	MW785 w/10ml	Tube: Neutralizing Buffer, 4mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83010 PDB NB	–	–	2583010PDBNB	–	

## ESK® SAMPLING KIT: BUFFERED PEPTONE WATER

Puritan	RELATED PRODUCTS				PRODUCT DESCRIPTION
	Copan	3M	Hardy Diagnostics	Medical Wire	
25-83004 PD BPW	–	–	2583004PDBPW	MW777 w/5ml	Tube: Buffered Peptone Water, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83004 PDB BPW	–	–	2583004PDBBPW	–	
25-83010 PD BPW	R4110 w/10ml	RS96010BPW w/10ml	2583010PDBPW	MW776 w/10ml	Tube: Buffered Peptone Water, 4 mL or 10 mL; Swab: Sterile polyester tip on polystyrene handle.
25-83010 PDB BPW	–	–	2583010PDBBPW	–	



CALL 800-321-2313 (USA & Canada) • 207-876-3311  
 EMAIL [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)  
 VISIT [puritanmedproducts.com](http://puritanmedproducts.com)

PURITAN MEDICAL PRODUCTS CO.  
 GUILFORD, MAINE, 04443, USA  
 ISO 13485:2016 CE



MADE IN THE USA